


讓你得到爆量的 **Plasmid**
以及肥滋滋的 **Colony**!



Terrific Broth

E. coli media for highest yields of plasmid DNA

Terrific broth (TB) is a powder formulation for increase bacterial growth and plasmid DNA.

1. 增加 4-7 倍 **Plasmid DNA**: Terrific broth 含有獨特的配方，跟 LB Broth 比起來，養的菌可以產生 **4-7 倍** 的 **plasmid**
2. 穩定的 pH 值，可以高濃度養菌：安心無負擔!
3. 保證：不含 **BSE** 成分
4. 特色：菌菌會很開心，養 **colony** 會長很大!



pUC19 plasmid with a 10kb insert:
TB (lane 1 and 2)
LB Broth (lane 3 and 4)

Cat. #	Product	Price	Cat. #	Product	Price
T20-114	Terrific Broth (TB), 500 g	2,500	fas-am-b	Fast-Media® Amp LB, 30 pouches	
GT-A01	Agar (for Bacteriology), 500 g	1,500	fas-am-s	Fast-Media® Amp Agar, 30 pouches	
A01-102	Agar (Molecular Biology Grade), 500 g	2,500	fas-am-x	Fast-Media® Amp XGal Agar, 20 pouches	
A1705	Agarose LE, 500 g		fas-kn-x	Fast-Media® kan XGal Agar, 20 pouches	
A2501	EZ Pack™ Agarose Tablets, 100 g	3,500	L12-111	LB Agar (MILLER'S), 500 g	3,000
GB0232	Glycerol, 500 ml		L12-112	LB Broth (MILLER'S), 500 g	1,500
AB0028	Ampicillin, sodium salt, 25 g		C03-104	Peptone, 500 g	2,500
KB0286	Kanamycin Sulfate, 25 g		T20-115	Tryptone, 500 g	2,500
TB0504	Tetracycline hydrochloride, 100 g		Y25-101	Yeast Extract, 500 g	2,500
ant-zn-1	Zeocin™ (solution), 1g		Y25-102	YPD Agar, 500g	
SL-90001	loop and needle, 500 pcs		Y25-103	YPD Broth, 500g	

F-TB20170615-1 版



萬造生技股份有限公司

Tel: (02)2908-3971 專線: 0988-753-268



Line ID: watson.biotechnology

e-mail: watson.biotechnology@gmail.com